
DIGITAL DISPLAY ALCOHOL BREATH TESTER

1. Breath inhaler
2. Power button
3. LCD display
4. Battery compartment
5. Key chain

Reference


BAC: Blood Alcohol Concentration
0.05% BAC: refers to there is 0.05g alcohol in each 100ml of blood

Features

1. Digital alcohol tester with audible alert
2. Range: 0.00-0.19 % BAC & 0.0-1.9 g/l
3. Quick response and resume
4. Auto power off
5. Use mouthpiece
6. Key chain
7. Use 2×“AAA” alkaline batteries

Installing battery

1. Slide open the battery compartment lid from the back of the unit.

2. Insert two “AAA” (UM-4) alkaline batteries with correct polarity as indicated.
3. Replace the battery compartment lid.


How to use the alcohol tester

Warm up

1. Press and hold the power button for 1 second, the alcohol tester LCD will turn on with one beep sound.

Precaution

- When after drinking, it is recommended to last after 20 minutes, it is because 20 minutes is approximately for alcohol to be absorbed into blood from the digestive organs, and residual alcohol remaining in the mouth takes this long to dissipate.
- For accurate result, do not repeat next testing within 3 minutes after last time.
- Do not smoke and eat for 20 minutes prior testing.
- Do not use chemicals to clean the tester directly.
- Do not blowing liquid into the breath pipe.
- Do not use any excitant substances near the tester, such as paint, insecticides, alcohol.
- Do not place the tester in the closed environment with contaminated air.
- Modifying or tampering with the unit's internal components may cause a malfunction.
- The unit is designed to be used as an aid to promote reasonable driving. It should not be used in anticipation of the use or operation of vehicles, machinery, or equipment. The unit makes no references whatsoever to blood alcohol content or other such measures of intoxication promulgated by government, medical or other agencies or bodies.


2. The “warm up” symbol is displayed, and the 10 to 00 countdown is started that indicate the tester is in warm-up mode.

Note


It always has a large number of gas & impurity adsorbed on the surface of alcohol tester, when first time taking out the tester from the packaging, or after long time no using for the tester, it need more time for cleaning and warm up the sensor, you should repress the power button and retry several times if possible.

Test

1. The “Blow” symbol is displayed when warm-up is succeed, and the 10 to 00 countdown is started that indicate the tester is in testing mode.

2. Blow into the breath inhaler for 3-5 seconds with normal power after a deep inhalation..
3. Read the test result from the alcohol tester LCD.
4. The “Caution” symbol is displayed if the alcohol concentration is the level:
 - ✓ 0.02% BAC (or 0.2g/l) to 0.05% BAC (or 0.5g/l)
 - 0.05% BAC (or 0.5g/l) to 0.08% BAC (or 0.8g/l)

5. The “Danger” symbol is displayed with beep audible alert if the alcohol concentration is the level:
 - ✓ Over & equal 0.05% BAC (or 0.5g/l)
 - Over & equal 0.08% BAC (or 0.5g/l)


6. If no blow into the breath inhaler at the test mode, the LCD will display “0.00 %BAC & 0.0g/l” value.

Note

To obtain more accurate result, manufacturer recommend you to check repeatedly 2 or 3 times.


Auto power off

After the result be display for about 10 seconds, the “OFF” symbol will appear for 2-3 seconds and the alcohol tester be turned off.

Replace mouthpiece

1. Contrarotate once, the mouthpiece will be loosened.
2. Replace the new mouthpiece on the right place, clockwise rotation, it will be ok.

Backlight

When the tester is turned on, the LCD backlight will light until tester auto power off.

Replace battery

The LCD flickering indicates low power of batteries, please replace two “AAA” alkaline batteries.

Specifications

- * Sensor: advanced semiconductor oxide alcohol sensor
- * Range: 0.00%-0.19% BAC (0.0-1.9g/l)
- * Accuracy: ±0.01% BAC (0.1g/l)
- * Power input: 2×1.5V “AAA” alkaline battery
- * Operating time: warm up: 10s, response: <5s
- * Auto power off: 1 minute approx
- * Dimension: 118×70×28mm
- * Environmental temperature: 5°Cæ°C

Notes

1. Regardless of whether the product is used, manufacturers and vendors will not bear any responsibility for the violations. For example: drunk driving and high-altitude operations.
2. Each person reacts differently to alcohol, so the test results are for reference only, which cannot serve as the criteria in determining the outcome.

Common Faults and Solutions

Faults	Possible Faults Reason	Solution
Power on failing	Put batteries incorrectly	Installing batteries in correct polarity
	Very low power batteries	Replacing the new alkaline batteries
	Circuit Fault	Contacting after-sale service
Blowing no response	Warm up not completing	Warm up for more times
	Circuit Fault	Contacting after-sale service
Blowing no result	Residual gas in tester	Warm up again
	Not blowing in schedule time	Blowing during “BLOW” countdown time
No blowing with test result	Residual gas in tester	Warm up for more times
Different result after several testing	Testing as soon as drinking	Testing 20 minutes after drinking
	Blowing time different	Deep breath for 3-5 seconds

Different country alcohol units

No.	BAC (Blood Alcohol Concentration)				BrAC (Breath Alcohol Concentration)	
	%	‰	g/l	mg/100ml	mg/l	ppm
1	0.01	0.1	0.1	10	0.05	25
2	0.02	0.2	0.2	20	0.10	50
3	0.03	0.3	0.3	30	0.15	75
4	0.04	0.4	0.4	40	0.20	100
5	0.05	0.5	0.5	50	0.25	125
6	0.06	0.6	0.6	60	0.30	150
7	0.07	0.7	0.7	70	0.35	175
8	0.08	0.8	0.8	80	0.40	200
9	0.09	0.9	0.9	90	0.45	225
10	0.10	1.0	1.0	100	0.50	250
11	0.11	1.1	1.1	110	0.55	275
12	0.12	1.2	1.2	120	0.60	300
13	0.13	1.3	1.3	130	0.65	325
14	0.14	1.4	1.4	140	0.70	350
15	0.15	1.5	1.5	150	0.75	375
16	0.16	1.6	1.6	160	0.80	400
17	0.17	1.7	1.7	170	0.85	425
18	0.18	1.8	1.8	180	0.90	450
19	0.19	1.9	1.9	190	0.95	475
20	0.20	2.0	2.0	200	1.00	500

Common units: %BAC-----USA
g/l-----European
mg/L-----Janpan
mg/100ml-----China

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